




Guardian® II Hemostasis Valve

Maintains hemostasis during interventions

The Guardian® II hemostasis valve is engineered with an 8F inner lumen and star-shaped seal technology to independently secure multiple devices for reduced blood loss during complex interventions.

Ergonomic Design


Promotes easy, single-handed operation

Rotolock

Provides option to secure devices


Star-Shaped Seal Technology

Designed to minimize blood loss while allowing independent guidewire movement




Guardian® II Hemostasis Valve Star-Shaped Seal

COPILOT® Valve Seal*



Available in Two Versions:



Ordering Information

Guardian® II Hemostasis Valve	
PRODUCT NO.	DESCRIPTION
8210	Guardian® II with Guidewire Insertion Tool
8211	Guardian® II with Guidewire Insertion Tool and VSI Torque Device
8215	Guardian® II NC with Guidewire Insertion Tool
8216	Guardian® II NC with Guidewire Insertion Tool and VSI Torque Device

Packaged in quantities of 25 units per box. Individual units are not sold separately. The Guardian® Hemostasis Valves are intended to maintain hemostasis during the use of diagnostic or interventional devices. Guardian® Hemostasis Valves are indicated for maintaining a seal around diagnostic/interventional devices with outside diameters up to 8F (0.105" or 2.67mm) during diagnostic/interventional procedures.

*Data on file at Teleflex. Comparative data may not be indicative of clinical performance.

1. Data on file at Teleflex. Comparative data may not be indicative of clinical performance.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüschi®, and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Poland +48 22 4624032

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, Arrow, Deknatel, Guardian, Hudson RCI, LMA, Pilling, Rüschi and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks and registered trademarks are property of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 09/2019.

© 2019 Teleflex Incorporated. All rights reserved. MC-005917 Rev 0