




Wattson™ Temporary Pacing Guidewire

Engineered for Efficiency

Changing the Pace of TAVR Procedures

The Wattson™ Temporary Pacing Guidewire offers dual functionality—a pigtailed guidewire that supports both valve delivery and ventricular bipolar pacing during transcatheter aortic valve replacement (TAVR) or balloon aortic valvuloplasty (BAV) procedures. The device offers a procedural alternative designed to help avoid a range of complications associated with traditional right ventricular pacing. From its flexible distal pigtail shape to its multiple electrode, bipolar design, the Wattson™ Temporary Pacing Guidewire is engineered to help reduce the risk of ventricular perforation while providing confidence in capture during rapid pacing.



2 in 1


TAVR guidewire +
pacing capabilities



Intraventricular
bipolar pacing

1st

First commercially available
bipolar temporary pacing
guidewire designed specifically
for TAVR and BAV



Ordering Information

Wattson™ Temporary Pacing Guidewire

MODEL	LENGTH	SHAFT O.D.	DISTAL PIGTAIL O.D.	GUIDEWIRE COATING
2250	280 cm	0.035" (0.89 mm)	30 mm	Lubricious Hydrophilic Coating (Distal 80cm)

The Wattson™ Temporary Pacing Guidewire is intended to introduce and position catheters and other interventional devices within the chambers of the heart, including those used within transcatheter aortic valve replacement (TAVR) procedures and balloon aortic valvuloplasty (BAV), while transmitting an electrical signal from an external pulse generator to the heart. The temporary pacing guidewire is not intended to remain in place following the clinical procedure.

The Wattson™ Temporary Pacing Guidewire is contraindicated for use in the coronary arteries and in the cerebrovasculature.

1. Bench test data completed by and on file at Teleflex. Bench test results may not necessarily be indicative of clinical performance.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation—a relentless pursuit of identifying unmet clinical needs—to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular access, interventional cardiology and radiology, anesthesia, emergency medicine, surgical, urology and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow™, Deknatel™, LMA™, Pilling™, QuikClot™, Rüschi™, UroLift™, and Weck™—trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Poland +48 22 4624032

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, Arrow, Deknatel, LMA, Pilling, QuikClot, Rusch, UroLift, Wattson, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. Refer to the Instructions for Use for a complete listing of the indications, contraindications, warnings, and precautions. Not all products may be available in all countries. Please contact your local representative.

© 2024 Teleflex Incorporated. All rights reserved. MC-006169 Rev 1. Revised 04/2024.